



member of SWEDENCARE

INNOVET ITALIA SRL  
Via Leonardo da Vinci 3  
I-35030 SACCOLONGO (PD)  
Tel. 0039 049 80 15 583  
Email [innovet@innovet.it](mailto:innovet@innovet.it)  
Web [www.innovet.it](http://www.innovet.it)

DATA SHEET

## Repy<sup>®</sup> Dress HA

Last update: 19/06/2026

### Sterile single-use dressing to promote the healing of skin lesions.

Thanks to the presence of Manuka honey and Adelmidrol, Repy<sup>®</sup> Dress HA promotes proper wound healing of lesions of any origin (e.g. traumatic, surgical, burns, avulsions, scratching).

Repy<sup>®</sup> Dress HA is a device for veterinary use.

Species: **Dog • Cat**



### Packaging and Prices

Repy<sup>®</sup> Dress HA 10x10 cm

10 units

70,50 €

cod. 2183



### Use

Clean the affected area and apply a thin layer of Repy<sup>®</sup> Gel before applying the dressing.

Open the package and remove one of the two protective films.

Apply the side of the gauze without the film to the affected skin surface, then remove the second protective film.

Secure the dressing with a self-adherent bandage.

Replace daily in case of heavy exudate; reduce the frequency of replacement as the re-epithelialisation process progresses and the exudates decrease. Replace at least every 3 days (or sooner if the dressing becomes completely dry).

Apply a thin layer of Repy<sup>®</sup> Gel each time the dressing is replaced.

Remove the dressing slowly and gently, using tweezers or two fingers (preferably while wearing sterile gloves). If the dressing sticks to the area, moisten it with sterile saline solution before removal.

### WARNINGS

Do not use if the packing has been opened or damaged, since the dressing is no longer considered sterile. Keep the product in a cool, dry place at room temperature.

### Properties

Repy<sup>®</sup> Dress HA is a sterile single-use gauze made of 100% pure cotton, impregnated with a solution of medical-grade Manuka honey and Adelmidrol. Monofloral Manuka honey helps maintain the proper moisture level during the healing process. Additionally, it controls the growth of surface microflora due to its high content of methylglyoxal (MGO), a natural antimicrobial agent in honey that has proven effective in reducing bacterial biofilm. Adelmidrol, a water-soluble analogue of palmitoylethanolamide, optimizes wound healing and improves the elasticity of the skin. It also provides a soothing and anti-redness effect.

### Functional principle

Natural cotton gauze

Manuka honey (98%\*)

Adelmidrol (2%\*)

\* of the impregnating solution

### Ingredients

Manuka honey (98%), Adelmidrol (2%)