



member of SWEDENCARE

INNOVET ITALIA SRL
Via Leonardo da Vinci 3
I-35030 SACCOLONGO (PD)
Tel. 0039 049 80 15 583
Email innovet@innovet.it
Web www.innovet.it

DATA SHEET

Repy® Gel

Last update: 02/05/2026

Re-epithelialisation promoting gel

Repy® Gel can be used to improve the functional and aesthetic outcome of skin re-epithelialisation in dogs and cats.

Repy® Gel isn't a veterinary drug. It's a dermatological product for dogs and cats.

Species: **Dog • Cat**



Packaging and Prices

Repy® Gel	75 ml	21,50 €	cod. 2101	
Repy® Gel	30 ml	15,30 €	cod. 2103	

Use

Apply a proper amount of product on the affected skin. Repeat two or more times daily.

Properties

Adelmidrol - a congener of palmitoylethanolamide - optimizes re-epithelialisation and improves skin elasticity. It also exerts a soothing and anti-redness effect. **Hyaluronic** and **traumatic acid** maintain integrity, hydration and elasticity of the skin. The **complex of plants extracts** (Echinacea purpurea and Usnea barbata) and **snail secretion** help to maintain the balance of the skin surface microflora. The addition of **Neem Oil** gives the product a repellent action against licking and insects (e.g., flies).

Functional principle

Aldemidrol

Hyaluronic acid

Traumatic acid

Echinacea

Usnea barbata

Snail Secretion Filtrate

Neem Oil

Ingredients

Aqua, echinacea purpurea extract, propylene glycol, melia azadirachta seed oil, snail secretion filtrate, algin, adelmidrol (INN), peg-40 hydrogenated castor oil, phenoxyethanol, ethylhexylglycerin, usnic acid, magnesium hyaluronate trans-2-dodecenedioate, sodium hydroxide.